

संदर्भ : कें.मु.वि-/जी-73

23 जुलाई 2018

**विषय:** IS 2062, IS 10748, IS 2879, IS 3975, IS 11513, IS 7466 के प्रॉडक्ट मानुयल

सक्षम अधिकारी द्वारा अनमोदित, निम्नलिखित प्रॉडक्ट मानुयल भा मा ब्यूरो वैबसाइट पर "Product Specific Information" के अंतर्गत होस्ट किए गए हैं:

([http://www.bis.org.in/cert/prod\\_specific\\_info.asp](http://www.bis.org.in/cert/prod_specific_info.asp))

- i. IS 2062:2011 - Hot Rolled Medium and High Tensile Structural Steel
- ii. IS 10748:2004 - Hot rolled steel strip for welded tubes and pipes
- iii. IS 2879:1998 - Mild Steel for Metal Arc Welding Electrodes
- iv. IS 3975:1999 - Low Carbon Galvanized Steel Wires, Formed Wires And Tapes For Armouring Of Cables
- v. IS 11513:2017 - Hot Rolled Carbon Steel Strip for Cold Rolling Purposes
- vi. IS 7466 : 1994 - Rubber Gaskets for Pressure Cookers

सभी संबन्धित से अनुरोध है कि इस का कार्यान्वयन सुनिश्चित करें।

(आदित्यदास)  
वैज्ञानिक सी (कें.मु.वि.-2)

प्रमुख, (कें.मु.वि.-2)

सभी क्षेत्रीय उप्महानिदेशकों/क्षेत्रीय / शाखा कार्यालयों को intranet माध्यम से परिचालित

प्रतिलिपि: ITS - इंटरनेट पर अपलोड करने के लिए

#### CENTRAL MARKS DEPARTMENT -2

Our Ref: CMD-2/G-73

23 July 2018

**Subject: Product Manuals for IS 2062, IS 10748, IS 2879, IS 3975, IS 11513, IS 7466**

Product Manuals, duly approved by the Competent Authority, have been hosted on BIS website under "Product Specific Information" ([http://www.bis.org.in/cert/prod\\_specific\\_info.asp](http://www.bis.org.in/cert/prod_specific_info.asp))

- i. IS 2062:2011 - Hot Rolled Medium and High Tensile Structural Steel
- ii. IS 10748:2004 - Hot rolled steel strip for welded tubes and pipes
- iii. IS 2879:1998 - Mild Steel for Metal Arc Welding Electrodes
- iv. IS 3975:1999 - Low Carbon Galvanized Steel Wires, Formed Wires And Tapes For Armouring Of Cables
- v. IS 11513:2017 - Hot Rolled Carbon Steel Strip for Cold Rolling Purposes
- vi. IS 7466 : 1994 - Rubber Gaskets for Pressure Cookers

All concerned are requested to ensure implementation of the same.

(Aditya Das)  
Sc. C (CMD-2)

Head, CMD-2

Circulated to DDGRs/ROs/BOs through intranet

Copy to: ITS - for hosting on Intranet